

How to Backup




identiMetrics™

Cheat Sheet


It's easy and important to back up your finger scan database. These instructions show you how to back up.

To backup your database, click on the Server Management Dropdown at the top of the Enroll Manager Window and select Backup.



This function creates a backup of the database and places it in the Requested Backup Folder on the server.

After selecting this function the Backup Dialog is displayed. You can use the Backup Folder name provided or you can enter you own.



Click on OK to run the backup or click on Cancel to close the dialog and quit the Backup function.

Need assistance? It's easy to get in touch with one of our Biometric Experts!



GO  .com



support@identimetrics.net



215-836-5640